

Laserglow Product Datasheet

LQS-1313 Passively Q-Switched Laser System

Laserglow Part Number: QD20405SX

This model is listed as **inactive** in our product database. Stock may be limited, and availability is subject to change without notice.

Similar Products:

For information about the other lasers in this product family visit:

http://www.laserglow.com/QD2

Ordering:

Order Online Now or Request Quote:

http://www.laserglow.com/QD20405SX

Series Specifications:

| Nominal Wavelength | 1313 nm |
|--------------------|------------|
| Output Type | Q-Switched |
| Laser Source Type | DPSS |



Overview:

The LQS-1313 Series of Diode-Pumped Solid-State (DPSS) Q-Switched Lasers are ideal for applications requiring a very high peak power or short pulse duration at 1313 nm.

These lasers are commonly used for Raman spectroscopy, material processing, and a broad range of other applications. The driver is available as a plugand-play benchtop system or an O.E.M. component designed for system integration.

Key Features:

- Pulse energy of 5 uJ 40 uJ
- Pulse repetition rate of 1 Hz 5 kHz
- Pulse duration in the range of 10-25 ns
- · Air cooled
- Runs on standard AC power (85 264 V, 47 63 Hz)
- 10,000 hour maintenance-free operating life (Expected)
- FDA/CDRH compliant Class IV enclosure

Package Includes:

- Laser Head
- · Driver/Power Supply

- Power Cable
- BNC Connector (LabSpec models only)
- Keys, Safety Interlock
- Hard-shell Carrying Case

Specifications:

This spec sheet has been generated specifically for part number QD20405SX, per your request, and data for the entire series is also displayed for your reference. The specs which are specific to QD20405SX have been highlighted below in **red + bold**.

| Output Power (mW) | >25 | >50, >100, >200 | |
|--|-------------------------------|---------------------------|--|
| Single Pulse Energy (µJ) | 5 | 10, 20, 40 | |
| Optimal Pulse Frequency (Hz) | 5000 | 5000 | |
| Output Power Stability (%RMS/4h) | <1, <3, <5 | <1, <3, <5 | |
| Wavelength Tolerance (+/- nm) | 1 | 1 | |
| Divergence (mrad, full angle) | <1.5 | <2 | |
| Beam Dimensions (mm, 1/e²) | 1.5 | 3 | |
| Warm-up Time (minutes) | 10 | 10 | |
| Avg. Pulse Duration (ns) | 10 | 25 | |
| Approximate Peak Power (W) | 500 | 1000, 4000 | |
| M² | <1.2 | <1.5 | |
| Operating Temperature Range (°C) | 10 to 35 | 10 to 35 | |
| Max. TTL Modulation Freq. (Hz) | 20000 | 20000 | |
| Minimum Pulsing Frequency (Hz) | 1 | 1 | |
| Modulation Input Signal | 0-5 VDC | 0-5 VDC | |
| Max. Power Input Duty Cycle | 100% | 100% | |
| Standard Warranty (months) | 12 | 12 | |
| MTTF (operational hours) | 10000 | 10000 | |
| Weight of Product or Laser Head (kg) | 0.6 | 0.9 | |
| Beam Height from Base Plate (mm) | 24.8 | 29 | |
| Dimensions of Product or Laser Head (mm) | 140.8 (l) x 73 (w) x 46.2 (h) | 155 (I) x 77 (w) x 60 (h) | |

CW: All specifications are based on performance at full output power and after the specified warmup period. Output characteristics may change if the laser is run at a different power level.

Q-Switched: Specifications are based on the laser pulsing at the specified design frequency. Output characteristics may change if the laser is run at a different frequency.

Power Supply Options:

These lasers are available with several different power supply options. The model that you have selected is highlighted below, and any other options are shown for easy reference.

| | Power Supply Type: | SM | SH |
|------------------------|-----------------------------|-------------------------------|--------------------------------|
| FDA-Compliant Standard | Input Power | 85v to 264v | 85v to 264v |
| | Power Supply Weight (kg) | 1.2 | 2.3 |
| | Dimensions (mm) | 133 (l) x 130 (w) x 65 (h) | 238 (l) x 146 (w) x 102 (h) |

^{*}Power supply may not be exactly as shown, see dimensional drawings on next 2 pages.


Regulatory Classification:

The model you have selected (QD20405SX) requires the following safety label(s):




^{*}Dimensions for fiber-integrated (I_) include laser head packaged inside.

Dimensional Drawing - Laser Form Factor: H:



Dimensional Drawing - Power Supply Form Factor: SH:



Accessories:

The most popular accessories for model QD20405SX are shown below. For additional details regarding these or other accessories please see our website or contact us directly.

| Part Number | Description | |
|-------------|---|---------------------|
| ACFLIRHXA | FC/PC Fiber Coupler/Collimator for IR wavelengths (1300 to 1600 nm) (installed and aligned) 11 mm diameter input lens Full Details: www.laserglow.com/ACF | |
| ACSLIRHXA | SMA-905 Fiber Coupler/Collimator for IR wavelengths (1300 to 1600 nm) (installed and aligned) 11 mm diameter input lens Full Details: www.laserglow.com/ACS | |
| ACALBMXXX | Carrying Case-102 Holds Lab/OEM M, R and S size, standard or LabSpec laser Full Details: www.laserglow.com/ACA | Included With Laser |
| AFF2002XX | Armored Fiber With FC/PC Connectors 200um Core Multimode 2m length Full Details: www.laserglow.com/AFF | |
| AFS2002XX | Armored Fiber With SMA 905 Connectors 200um Core Multimode 2 m length Full Details: www.laserglow.com/AFS | |
| AGFMIR4XX | LSG-MIR-NF-4 Fit-Over Safety Goggles Mid-IR Range Output: OD 4+ at 945-10600 nm CE Certified Full Details: www.laserglow.com/AGF | |

FOR MORE INFORMATION PLEASE CONTACT:

LASERGLOW TECHNOLOGIES
99 Ingram Dr. Unit B, North York, ON, Canada M6M2L7
Tel. (416) 729-7976 Fax (716) 322-3510
sales@laserglow.com www.laserglow.com

E&OE: Data included in this sheet may be subject to change without notice.

Please confirm critical specifications with our staff prior to ordering.